


Compulink SPEX UPC CD Installation Guide

Insert SPEX UPC CD

1. After putting the SPEX UPC disc in your computer, **open** “Computer” or “My Computer” to see what drive the disc is on. **Keep note of the drive letter**, you’ll need it later!




Loading Data to Compulink

1. **Open** Compulink.


2. Click **Utility -> Miscellaneous Database Utilities -> Optical Utilities -> Frames Data Update**



3. Remember the drive letter from Step 1? Enter the letter and **click "OK"**.



4. **Click** “OK” in the pop-up window that appears.



5. The data should now be loading into Compulink. This might take about 20 minutes or so to complete.



Congratulations! You have now successfully loaded Frames Data SPEX UPC CD into Compulink.

If you have any questions please contact technical support at techsupport@framesdata.com or by phone at 1-800- 821-6069 ext.1.